

10/535442

JC06 Rec'd PCT/PTO 19 MAY 2005

WO 2004/046379

PCT/FI2003/000888

SEQUENCE LISTING

<110> MoBiDiag Oy

<120> Nucleic acid probes, broad-range primers, and methods in which they are used

<130> 2022264FI

<160> 77

<170> PatentIn version 3.1

<210> 1

<211> 23

<212> DNA

<213> Streptococcus pneumoniae

<400> 1  
ctcaaaagaa ggtcttcacc atc 23

<210> 2

<211> 23

<212> DNA

<213> Streptococcus pneumoniae

<400> 2  
gccatattca agttttttatt gag 23

<210> 3

<211> 23

<212> DNA

<213> Streptococcus pneumoniae

<400> 3  
agccagatga ttcgattact gtt 23

<210> 4

<211> 23

<212> DNA

<213> Streptococcus pneumoniae

<400> 4  
ttgagaccgt ctttacagtc ctt 23

<210> 5

<211> 23

<212> DNA

<213> Streptococcus pyogenes

<400> 5  
gttttgctc tcatattaaa gtc 23

<210> 6  
<211> 23  
<212> DNA  
<213> Streptococcus pyogenes

<400> 6  
tcctttattga agcagataat tcc 23

<210> 7  
<211> 23  
<212> DNA  
<213> Streptococcus pyogenes

<400> 7  
cgttgaaaca gtttttacag tct 23

<210> 8  
<211> 23  
<212> DNA  
<213> Chlamydia pneumoniae

<400> 8  
ggtcttcatc atctagtcta tga 23

<210> 9  
<211> 23  
<212> DNA  
<213> Chlamydia pneumoniae

<400> 9  
ggttgtagac wacagcattg acg 23

<210> 10  
<211> 23  
<212> DNA  
<213> Chlamydia pneumoniae

<400> 10  
agccatggca gttattgct cta 23

<210> 11  
<211> 23  
<212> DNA  
<213> Chlamydia pneumoniae

<400> 11  
ggattgatgt tsgcatttta gag 23

<210> 12  
<211> 23  
<212> DNA  
<213> Chlamydia pneumoniae

<400> 12  
gggtattgtc wtcgtagata atg 23

<210> 13  
<211> 22  
<212> DNA  
<213> Chlamydia pneumoniae

<400> 13  
ttattgctct aggattgatg tt 22

<210> 14  
<211> 23  
<212> DNA  
<213> Chlamydia pneumoniae

<400> 14  
gcattttaga gsacgggggt att 23

<210> 15  
<211> 24  
<212> DNA  
<213> Mycoplasma pneumoniae

<400> 15  
tcaaactaac ccttaaagac aact 24

<210> 16  
<211> 23  
<212> DNA  
<213> Mycoplasma pneumoniae

<400> 16  
tgaaacagtg tttacggtac tcc 23

<210> 17  
<211> 23  
<212> DNA  
<213> Haemophilus influenzae

<400> 17  
atgataattc tgtatcggtg caa 23

<210> 18  
<211> 23  
<212> DNA  
<213> Haemophilus influenzae

<400> 18  
agcagaagtt attatgactg tgc 23

<210> 19  
<211> 23  
<212> DNA  
<213> Haemophilus influenzae

<400> 19  
gacgataact cttataaaagt atc 23

<210> 20  
<211> 24  
<212> DNA  
<213> Haemophilus influenzae

<400> 20  
gataccgatg acggtactgg tttg 24

<210> 21  
<211> 23  
<212> DNA  
<213> Staphylococcus aureus

<400> 21  
agtgtgggaa attgtcgata ata 23

<210> 22  
<211> 25  
<212> DNA  
<213> Staphylococcus aureus

<400> 22  
tgaagttggtt attgaaaaag ataac 25

<210> 23  
<211> 23  
<212> DNA  
<213> Staphylococcus aureus

<400> 23  
ggattaaagt aacggataac gga 23

<210> 24  
<211> 24  
<212> DNA  
<213> Staphylococcus aureus

<400> 24  
ctgtcgaagt tatttttaact gttt 24

<210> 25  
<211> 23  
<212> DNA

<213> Staphylococcus aureus

<400> 25

cgacttcaga gagaggtttg cac

23

<210> 26

<211> 23

<212> DNA

<213> Pseudomonas aeruginosa

<400> 26

gaaatcagca tcaccatcca tac

23

<210> 27

<211> 23

<212> DNA

<213> Pseudomonas aeruginosa

<400> 27

atacggatga gtcgatcact gtc

23

<210> 28

<211> 23

<212> DNA

<213> Pseudomonas aeruginosa

<400> 28

gacgacaaca cctacaaggt gtc

23

<210> 29

<211> 24

<212> DNA

<213> Pseudomonas aeruginosa

<400> 29

gacaccgacg atggcaccgg tctg

24

<210> 30

<211> 23

<212> DNA

<213> Neisseria gonorrhoeae

<400> 30

ttcgacaaca acagctacaa aat

23

<210> 31

<211> 23

<212> DNA

<213> Neisseria gonorrhoeae

<400> 31

atatggtgtt tgaagtattg gac

23

<210> 32  
<211> 23  
<212> DNA  
<213> *Neisseria gonorrhoeae*

<400> 32  
gacaaaatca cggtaacgat aca 23

<210> 33  
<211> 23  
<212> DNA  
<213> *Neisseria gonorrhoeae*

<400> 33  
gacaacaaca gctacaaaat ctc 23

<210> 34  
<211> 23  
<212> DNA  
<213> *Escherichia coli*

<400> 34  
tattcgaggt ggtagataac gct 23

<210> 35  
<211> 23  
<212> DNA  
<213> *Escherichia coli*

<400> 35  
ttcacgccga taactctgtc tct 23

<210> 36  
<211> 23  
<212> DNA  
<213> *Escherichia coli*

<400> 36  
cgccgataac tctgtctctg tac 23

<210> 37  
<211> 23  
<212> DNA  
<213> *Escherichia coli*

<400> 37  
gacgataact cctataaagt gtc 23

<210> 38  
<211> 24

<212> DNA  
<213> Escherichia coli

<400> 38  
gacacggatg acggcaccgg tctg

24

<210> 39  
<211> 23  
<212> DNA  
<213> Moraxella catarrhalis

<400> 39  
agctgccgag gttattatga cgg

23

<210> 40  
<211> 23  
<212> DNA  
<213> Moraxella catarrhalis

<400> 40  
gatgataatt catacaaagt atc

23

<210> 41  
<211> 23  
<212> DNA  
<213> Moraxella catarrhalis

<400> 41  
tgtggatatc caccctgaag aag

23

<210> 42  
<211> 23  
<212> DNA  
<213> Moraxella catarrhalis

<400> 42  
ataccgatga tggtacaggc ttg

23

<210> 43  
<211> 23  
<212> DNA  
<213> Legionella pneumophila

<400> 43  
gatgataatt cctacaaggt atc

23

<210> 44  
<211> 23  
<212> DNA  
<213> Legionella pneumophila

<400> 44

agcagccgaa gtcacatga cag 23

<210> 45  
<211> 23  
<212> DNA  
<213> Legionella pneumophila

<400> 45  
tgtagatatt cataaagaag aag 23

<210> 46  
<211> 23  
<212> DNA  
<213> Legionella pneumophila

<400> 46  
atacagatga tggaaccggt ttg 23

<210> 47  
<211> 23  
<212> DNA  
<213> Legionella pneumophila

<400> 47  
cagatgatgg aaccggtttg cat 23

<210> 48  
<211> 23  
<212> DNA  
<213> Fusobacterium necrophorum

<400> 48  
caacatcagc caggggatta cat 23

<210> 49  
<211> 23  
<212> DNA  
<213> Fusobacterium necrophorum

<400> 49  
ggaattccag tagacatata tcc 23

<210> 50  
<211> 24  
<212> DNA  
<213> Fusobacterium necrophorum

<400> 50  
tgaaaatgat aactataaag tgtc 24

<210> 51



<211> 23  
<212> DNA  
<213> Staphylococcus aureus

<400> 51  
cggtaacgaa atagatgtaa caa 23

<210> 52  
<211> 23  
<212> DNA  
<213> Staphylococcus aureus

<400> 52  
agaagataat ggacgtggta tgc 23

<210> 53  
<211> 23  
<212> DNA  
<213> Staphylococcus aureus

<400> 53  
gtaaaccgac agtcgaagtt atc 23

<210> 54  
<211> 24  
<212> DNA  
<213> Staphylococcus aureus

<400> 54  
caactgataa acggggatta catc 24

<210> 55  
<211> 23  
<212> DNA  
<213> Staphylococcus aureus

<400> 55  
ggacaaggcg gctataaaac ttc 23

<210> 56  
<211> 23  
<212> DNA  
<213> Streptococcus pyogenes

<400> 56  
tggagatgat attaaggttg tta 23

<210> 57  
<211> 23  
<212> DNA  
<213> Streptococcus pyogenes

<400> 57  
gtcagtgtgg cagatagcgg acg 23

<210> 58  
<211> 23  
<212> DNA  
<213> Streptococcus pyogenes

<400> 58  
ggattccac cgttcaagtt att 23

<210> 59  
<211> 23  
<212> DNA  
<213> Streptococcus pyogenes

<400> 59  
caacagatgc tacgggattg cac 23

<210> 60  
<211> 23  
<212> DNA  
<213> Streptococcus pyogenes

<400> 60  
ggtcagggtg gctacaaaac gtc 23

<210> 61  
<211> 23  
<212> DNA  
<213> Streptococcus pneumoniae

<400> 61  
tggtgatcgt attgatgtaa cta 23

<210> 62  
<211> 23  
<212> DNA  
<213> Streptococcus pneumoniae

<400> 62  
ctaacggttc aagaccatgg acg 23

<210> 63  
<211> 23  
<212> DNA  
<213> Streptococcus pneumoniae

<400> 63  
gaattccaac tggtgaggtt atc 23

<210> 64  
<211> 23  
<212> DNA  
<213> Streptococcus pneumoniae

<400> 64  
cgaccgatgg cgctggtctt cat 23

<210> 65  
<211> 23  
<212> DNA  
<213> Streptococcus pneumoniae

<400> 65  
ggtcaagggtg gctataagac atc 23

<210> 66  
<211> 23  
<212> DNA  
<213> Mycoplasma pneumoniae

<400> 66  
tgctaatact attgccgttg ttt 23

<210> 67  
<211> 23  
<212> DNA  
<213> Mycoplasma pneumoniae

<400> 67  
ataactgtta gtgacaacgg tgc 23

<210> 68  
<211> 23  
<212> DNA  
<213> Mycoplasma pneumoniae

<400> 68  
agattttctac gattgacacc gtc 23

<210> 69  
<211> 23  
<212> DNA  
<213> Mycoplasma pneumoniae

<400> 69  
gataacgatt cttataagat tgc 23

<210> 70  
<211> 294  
<212> DNA  
<213> Moraxella catarrhalis

<400> 70  
 cgaccagggga tgtatattgg tgataccgat gatggtacag gcttgcacca tatggtgttt 60  
 gaggtggtgg ataatgccat tgatgaggca ttggcaggtc actgtgatga gattaatatt 120  
 atcgtccatg acgatgaatc tgtttcgggtg atggactatg ggcgtgggat tcctgtggat 180  
 atccaccctg aagaagggtg atcagctgcc gaggttatta tgacgggtgct tcatgcaggc 240  
 ggtaaatttg atgataattc atacaaagta tctggggggc tgcacggcgt agga 294

<210> 71  
 <211> 257  
 <212> DNA  
 <213> Legionella pneumophila

<400> 71  
 ggggatacag atgatggaac cggtttgcac cacatggttt ttgaggttgt agataattca 60  
 atagatgagt ctctggcagg atattgcaag gaaatTTTTTg ttaccatcca tagcgatgag 120  
 tcaattacag ttaaggacga tggccgtggg attcctgtag atattcataa agaagaaggc 180  
 aaatcagcag ccgaagtcac catgacagtc ctacatgctg gaggtaaatt tgatgataat 240  
 tcctacaagg tatctgg 257

<210> 72  
 <211> 290  
 <212> DNA  
 <213> Fusobacterium necrophorum

<400> 72  
 cgcccagggga tgtacatcgg aacaacatca gccaggggat tacatcattt agtatgggaa 60  
 gtggtagata attctgtgga tgaagcattt gctgggtatt gtaataggat tactgtaagt 120  
 attttgcttg ataacattat tcaagtagag gataatggaa gaggaattcc agtagacata 180  
 catccaaaat atggaaaatc cgctttggaa attgtattga cggattaca tgctggggga 240  
 aaatttgaaa atgataacta taaagtgtca ggtggactgc acggagttgg 290

<210> 73  
 <211> 291  
 <212> DNA  
 <213> Streptococcus mitis

<220>  
 <221> misc\_feature  
 <222> (18)..(18)  
 <223> n is a, t, c, or g

<220>  
 <221> misc\_feature  
 <222> (288)..(288)  
 <223> n is a, t, c, or g

<400> 73  
 agtccccggga tgtacatngg gtcaacttca aaagaaggtc ttcaccatct agtctgggaa 60  
 attggtgata actcaattga cgaggccttg gcaggatttg ctagccatat tcaagtcttt 120  
 attgagccag ataattcgat tacagttggt gatgatggac gtggtatccc agtcgatatt 180  
 caggaaaaaa caggtcgacc tgccggtgag actgtcttta ctgttcttca cgctggagga 240  
 aagtttggcg gtggcggata taaggctctca gacggctctc atggcgtngg a 291

<210> 74  
 <211> 288  
 <212> DNA  
 <213> Streptococcus oralis

<400> 74  
 agaccgggga tgtatatagg atcgaccgac ggtgctggtc tccatcatct agtctgggaa 60  
 atcgtggata atgcggttga cgaagccttg tctggatttg gtgatcgcat cgatgtgacg 120  
 attaataagg acgggagttt aacggttcaa gaccacggac gtgggatgcc aacgggaatg 180  
 cagcccatgg gaattccaac tggtgaagtt atctttacca ttctccacgc tggagggaaa 240  
 ttcggtcaag gtggctataa gacatctggt ggtctgcatg gggttgga 288

<210> 75  
 <211> 294  
 <212> DNA  
 <213> Haemophilus parainfluenzae

<400> 75  
 agaccgggga tgtatatcgg ggataccgat gatggaacag gcctacacca tatggtgttt 60  
 gaggtggtgg ataacgctat cgatgaagcg cttgctggct attgttccga tattatcgtc 120  
 actattcatg atgataattc tgtttccgta caagatgacg gccgcggtat tccggtagat 180  
 attcaccag aagaaggggt ttctgcggca gaagcaatca tgacagtact tcacgcaggt 240  
 ggtaaattcg atgataactc ttataaagta tcagggggac tacacggcgt tgga 294

<210> 76  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence

<220>

<223> primer

<220>

<221> misc\_feature

<222> (6)..(6)

<223> w i s a o r t

<220>

<221> misc\_feature

<222> (9)..(9)

<223> k i s g o r t

<220>

<221> misc\_feature

<222> (15)..(15)

<223> y i s c o r t

<220>

<221> misc\_feature

<222> (18)..(18)

<223> h i s a, c o r t

<400> 76

cgtccwggka tgtayathgg

20

<210> 77

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> primer

<220>

<221> misc\_feature

<222> (3)..(3)

<223> h i s a, c o r t

<220>

<221> misc\_feature

<222> (6)..(6)

<223> r i s a o r t

<220>

<221> misc\_feature

<222> (9)..(9)

<223> r i s a o r t

<220>

<221> misc\_feature

<222> (12)..(12)  
<223> w is a or t

<220>  
<221> misc\_feature  
<222> (15)..(15)  
<223> w is a or t

<220>  
<221> misc\_feature  
<222> (18)..(18)  
<223> d is a, g or t

<400> 77  
cchacrc crt gwaawccdcc

20